

Colorado Needle Microdissection




Specialties and procedures

Plastic Surgery

- Abdominoplasty
- Bi-coronal flap incision
- Blepharoplasty
- Breast Augmentation
- Cleft Lip/Cleft Palate
- Endoscopic brow lift
- Face Lift
- Otoplasty
- Platysma band notching
- Rhinoplasty

Neuro Surgery

- Coronal Incision
- Laminectomy
- Skull Base
- Spin/Rhizotomy



Specialties and procedures

Oral Maxillofacial Surgery

- Coronal Incision
- Facial incision trauma/plating
- Frenectomy
- Gingivectomy
- Intraoral Incision
- Palatal donor grafts
- Preauricular Incision
- Skin incision for hip graft
- Submandibular incision

Ear Nose and Throat Surgery

- ENT Cosmetic
- Neck Dissection
- Otoplasty
- Palatoplasty
- Preauricular Incision
- Rhinoplasty
- Tonsillectomies
- Turbinectomy
- Tympanomeatal Flap



In Service/Case Coverage/Trouble Shooting

The purpose of electrosurgery is to control bleeding₁

Two Modes:

Monopolar (Colorado Needle, single tip) - patient must have a grounding pad to complete circuit₃

- Used to access/open the surgical site

Bipolar (Forceps, 2x tipped) - no grounding pad needed because circuit completed between dual points₃

- Used once inside the surgical site

